

HCDOAV-154	5 Panel Cup AMP1000/COC300/OPI2000/PCP25/THC50 *CLIA WAIVED
HCDOAV-254	5 Panel Cup AMP1000/COC300/mAMP1000/OPI2000/THC *CLIA WAIVED ★
HCDOAF-274A6	7 Panel Cup w/ Adulteration Strip, AMP1000/BZO300/COC300/MET1000/OPI2000/OXY100/THC50 (OX/SG/pH) (GL/NI/CR) *CLIA WAIVED
HCDOAV-4104	10 Panel Cup AMP1000/BAR300/BZO300/COC300/MET1000/MTD300/OPI2000/OXY100/PCP25/THC50
IS2-DJA-1127-300	12 Panel Cup w/Adulteration Strip AMP500/BAR300/BZO300/COC150/MET500/MDMA500/MOR300/MTD300/OXY100/PCP25/TCA1000/THC50 (CRE, pH, OX)
HCODAV-6125A3A	12 Panel Cup AMP1000, BAR300, BUP10, BZO300, COC300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP25, THC50 (CR/SG/pH) *CLIA WAIVED
USSCUPA-12CLIA	12 Panel Cup w/Adulteration Strip AMP1000/BAR300/BZO300/COC300/MET1000/MDMA500/MOP300/MTD300/OXY100/PCP25/TCA1000/THC50 (CR,SG,PH) *CLIA WAIVED
US2-DJO-1137-300L	13 Panel Cup AMP500/BAR300/BUP10/BZO300/COC150/MET500/MDMA500/MOR300/MTD300/OXY100/PCP25/TCA1000/THC50
HCDOAV-8135EF1TA3A	13 Panel Cup, w/Adulteration Strip AMP1000/BUP10/BZO300/COC300/ETG500/FEN10/MDMA500/MET1000/MOP300/MTD300/OXY100/THC50/TRA200 (CR/SG/pH) (FUO)
HCDOAV-2145	14 Panel Cup, AMP1000/BAR300/BUP10/BZO300/COC300/EDDP300/MDMA500/MET1000/MOR300/MTD300/OXY100/PCP25/TCA1000/THC50

- Accurate
- Affordable
- Min. Vol. To Run is 5 mL
- Negative Results in > 1 minute
- Available in CLIA Waived 501K, Forensic Use, Insurance


★ - Best Seller



Multi-Drug Screen Cup

For Rapid Drug & Alcohol Screening

19525 Janacek Ct, Suite 104 Brookfield, WI 53045 ☎ +1 262.432.9308




- Accurate
- Affordable
- Min. Vol. To Run is 5 mL
- Negative Results in > 1 minute
- Available in CLIA Waived 501K, Forensic Use, Insurance

Directions for use


Step 1
Remove the cup from the foil pouch.




Step 2
Collect the urine in the cup.



The maximum amount of urine (120mL).





The minimum amount of urine required (5mL).



Step 3
Secure cap tightly.



Step 4
At 5 minutes, peel the label and read the results. Do not read after 60 minutes.



Negative(-)



Positive(+)



Invalid



Multi-Drug Screen Cup

For Rapid Drug & Alcohol Screening